

CG TAINER BIOTECH 250

Reusable multiway packaging
with a concept



SPECIFICATIONS & FACTS

Folding IBC	Tare weight: 32 kg	<ul style="list-style-type: none"> • The Bag in Box system: Shuttle container for sterile transport of biopharmaceutical process solutions • Transport safety confirmed with vibration test
variable filling volumes from 100 to 200 l	Hygienic surfaces	
Single-use bag*	Optimum use of space	
Dimensions: 80 cm x 60 cm x 83 cm	Extraction line in the bottom compartment	

Your challenge – our solution

DISMANTLING

In 4 steps

1 Loosen transport safety
(only for version < 200 l)




● Safety straps

2 Cover plate
(only for version < 200 l),
remove foam panel
and intermediate floor



● Red flap



● Cover plate

● Foam plate

● Intermediate floor

3 Dispose of single-use
bag after use*

4 Fold box



● CG Tainer

Please note: Return all components in one set up CG Tainer.


Up to 4 packages fit in one CG Tainer.

Other CG Tainers can simply be returned folded up.

* For safety and environmental reasons, we are not allowed to accept used single-use bags.



Optimally fastened



Safely sealed



Sterile packaged



Tamper-proof

WHAT CAN WE DO FOR YOU? WE ARE HERE FOR YOU.

CG 1962 | CG CHEMIAKLIEN GESELLSCHAFT MIT BESCHRÄNKTER HAFTUNG & CO. KG

🏠 Ulmerstraße 1, 30880 Laatzen

✉ info@cg-pharma.de

☎ +49 (0) 5118 7803-0

🌐 www.cg-pharma.de

CG

SEIT 1962

A company of the
CG GROUP

DISMANTLING

In 4 steps

RETURNS

What you return to us

1 Loosen transport safety
(only for version < 200 l)




● Safety straps

2 Cover plate
(only for version < 200 l),
remove foam panel
and intermediate floor



● Red flap



● Cover plate

● Foam plate

● Intermediate floor

3 Dispose of single-use
bag after use*

4 Fold box



● CG Tainer

Please note: Return all components in one set up CG Tainer.

Up to 4 packages fit in one CG Tainer.

Other CG Tainers can simply be returned folded up.

* For safety and environmental reasons, we are not allowed to accept used single-use bags.



Optimally fastened



Safely sealed



Sterile packaged



Tamper-proof

WHAT CAN WE DO FOR YOU? WE ARE HERE FOR YOU.

CG 1962 | CG CHEMIAKLIEN GESELLSCHAFT MIT BESCHRÄNKTER HAFTUNG & CO. KG

🏠 Ulmerstraße 1, 30880 Laatzen

✉ info@cg-pharma.de

☎ +49 (0) 5118 7803-0

🌐 www.cg-pharma.de

A company of the

